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Qwant

The search engine that respects your privacy.

ncbi blast

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 U.S. National Library of Medicine

 National Center for Biotechnology Information

BLAST[®]

HomeRecent ResultsSaved Str

Basic Local Alignment Search Tool

BLAST finds regions of similarity between biological sequences. The program compares nucleotide or protein sequences to sequence databases and calculates the statistical significance.

[Learn more](#)

NEWS


BLAST+ 2.10.1 is released – Fix for TBLASTN Multi-Threading issue.

This version supports pulling databases from our FTP site as well from [cloud providers](#) or our [BLAST+Docker solution](#).

Thu, 18 June 2020 12:00:00 EST

[More BLAST news...](#)


Web BLAST

Nucleotide BLAST
nucleotide ► nucleotide

blastx
translated nucleotide ► protein

tblastn
protein ► translated nucleotide

Protein BLAST
protein ► protein

blastn

blastp

blastx

tblastn

tblastx

BLASTN programs search nucleotide databases using a nucleotide query. [more...](#)

Enter Query Sequence

Enter accession number(s), gi(s), or FASTA sequence(s) [?](#)[Clear](#)Query subrange [?](#)

From

To

Or, upload file

Parcourir...

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Job Title

Echantillon inconnu n12

Enter a descriptive title for your BLAST search [?](#)☐ Align two or more sequences [?](#)

Choose Search Set

Database

☒ Standard databases (nr etc.): ☐ rRNA/ITS databases ☐ Genomic + transcript databases ☐ BetacoronavirusNucleotide collection (nr/nt) [?](#)

Organism

Optional

arthropods (taxid:6656)

☐ exclude[+](#)Enter organism common name, binomial, or tax id. Only 20 top taxa will be shown [?](#)

Exclude

Optional

☐ Models (XM/XP) ☐ Uncultured/environmental sample sequences

Limit to

Optional

☐ Sequences from type material

Entrez Query

Optional

 [YouTube](#) [Create custom database](#)Enter an Entrez query to limit search [?](#)

Program Selection

Optimize for

- ☒ Highly similar sequences (megablast)
☐ More dissimilar sequences (discontiguous megablast)
☐ Somewhat similar sequences (blastn)

Choose a BLAST algorithm [?](#)

BLAST

Search **database Nucleotide collection (nr/nt)** using **Megablast** (Optimize for highly similar sequences)☐ Show results in a new window

CTRL-V

?? idée

Clickkk !

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Job Title	Echantillon inconnu n12		
RID	M4DXJVG9014	Search expires on 08-24 17:32 pm	Download All ▼
Program	BLASTN ?	Citation	▼
Database	nt	See details	▼
Query ID	lcl Query_32393		
Description	Echantillon inconnu n12		
Molecule type	dna		
Query Length	420		
Other reports	Distance tree of results	MSA viewer	?

Filter Results

Organism only top 20 will appear ☐ exclude


[+ Add organism](#)

Percent Identity
 to

E value
 to

Query Coverage
 to

[Filter](#) [Reset](#)

- Descriptions
- Graphic Summary
- Alignments
- Taxonomy

Sequences producing significant alignments				Download	Manage Columns	Show	100	?
☒ select all 100 sequences selected				GenBank	Graphics	Distance tree of results		
	Description	Max Score	Total Score	Query Cover	E value	Per. Ident	Accession	
☒	Leptoneta convexa voucher 9040712 cytochrome c oxidase subunit I (COI) gene, partial cds; mitochondrial	776	776	100%	0.0	100.00%	JN378156.1	
☒	Leptoneta microphthalmia voucher 9040711 cytochrome c oxidase subunit I (COI) gene, partial cds; mitochondrial	593	593	100%	1e-166	92.14%	JN378155.1	
☒	Cybaeota nana cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial	407	407	99%	1e-110	84.21%	FJ263787.1	
☒	Erigone autumnalis haplotype 1 cytochrome c oxidase subunit 1 (COI) gene, partial cds; mitochondrial	405	405	99%	5e-110	84.29%	MK521227.1	
☒	Neoramia sp. CG178 isolate ARACG000178 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial	405	405	99%	5e-110	84.21%	KY017938.1	


Leptoneta convexa voucher 9040712 cytochrome c oxidase subunit I (COI) gene, partial cds; mitochondrial

Sequence ID: [JN378156.1](#) Length: **817** Number of Matches: **1**

Range 1: 1 to 420 [GenBank](#) [Graphics](#)

▼ [Next Match](#) ▲ [Previous Match](#)

Score	Expect	Identities	Gaps	Strand
776 bits(420)	0.0	420/420(100%)	0/420(0%)	Plus/Plus
Query 1	TCGTATGAATAAATTTGAGTTTTTGGTTATTACCCCCTTCTTTATTTTTATTAGTGGTTTC	60		
Sbjct 1	TCGTATGAATAAATTTGAGTTTTTGGTTATTACCCCCTTCTTTATTTTTATTAGTGGTTTC	60		
Query 61	TAGAATAGTTGAAATAGGAGTAGGAACAGGTTGGACTGTGTATCCTCCATTGTCTTCTAA	120		
Sbjct 61	TAGAATAGTTGAAATAGGAGTAGGAACAGGTTGGACTGTGTATCCTCCATTGTCTTCTAA	120		
Query 121	TGTAGGGCATGCTTTTGTGTCGGTTGATTTTGCTATTTTTTCTCTTCATTTAGCTGGGGC	180		
Sbjct 121	TGTAGGGCATGCTTTTGTGTCGGTTGATTTTGCTATTTTTTCTCTTCATTTAGCTGGGGC	180		
Query 181	TTCTTCTATTATAGGGGCTGTTAATTTTCATTACTACTATTATTAATATGCGAACTTTTGG	240		
Sbjct 181	TTCTTCTATTATAGGGGCTGTTAATTTTCATTACTACTATTATTAATATGCGAACTTTTGG	240		
Query 241	TATGAGTATGGAAAAGGTTCCCTTGTTTGTTTGATCAGTGTTAGTGACTGCTGTTTTGTT	300		
Sbjct 241	TATGAGTATGGAAAAGGTTCCCTTGTTTGTTTGATCAGTGTTAGTGACTGCTGTTTTGTT	300		
Query 301	GTTGTTATCTTTACCTGTATTAGCTGGTGCTATTACTATGTTGTAACTGATCGTAATTT	360		
Sbjct 301	GTTGTTATCTTTACCTGTATTAGCTGGTGCTATTACTATGTTGTAACTGATCGTAATTT	360		
Query 361	TAATACCTCttttttttGACCCTGCTGGCgggggggATCCAGTGTTATTTCAACATTTGTT	420		
Sbjct 361	TAATACCTCTTTTTTTGACCCTGCTGGCGGGGGGATCCAGTGTTATTTCAACATTTGTT	420		


Séquençage 'Sanger'

100-1200 bp


MinION 900 euro

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Consommables
539 euro

Ordi performant
3500 euro

G
A
C
T
G
A
A
G
C
T
G
T
T

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